

Bipolar Tonsillectomy with To-BiTE™ non-stick bipolar clamp



Bipolare Tonsillektomie mit der To-BiTE™ non-stick bipolaren Zange

To-BiTE™

FOUR FUNCTIONS

„Four functions in one instrument make tonsillectomies safer, faster, easier.“

P. Tolsdorff, MD, Bad Honnef (Germany)

1 Coag
Koagulieren

2 Dissect
Spreizen

3 Aspirate
Saugen

4 Grasp
Greifen



Made in Germany.



Tonsillectomies can be performed very quickly with little, but precise and effective bipolar coagulation while most of the dissection is “cold” or “blunt”. The glossopharyngeal nerve can be identified and spared.¹



Tonsillar Dissection.
Präparieren der Tonsille.



Wound immediately after tonsillectomy.
Wunde unmittelbar nach der Tonsillektomie.

¹ Tolsdorff P. Radiofrequency Micro-Bipolar Tonsillectomy. ENT & audiology news. 2010; Vol 18 No. 6: 14




+



=

Bipolar
Radiofrequency
Tonsillectomy



Ordering Information Bestellinformation

To-BiTE™ non-stick bipolar clamp - non-stick bipolare Zange

70 09 60 SG	To-BiTE™ non-stick bipolar clamp To-BiTE™ non-stick bipolare Zange
37 01 54 R	Bipolar (CURIS®) cable for To-BiTE™, length: 3 m Bipolares (CURIS®) Kabel für To-BiTE™, Kabellänge: 3 m

CURIS® 4 MHz Radiofrequency Generator - 4 MHz Radiofrequenz-Generator

87 00 10 CURIS® basic set with single-use patient plates | Basis-Set mit Einwegneutralelektrode:

36 01 00-01	CURIS® 4 MHz radiofrequency generator (incl. mains cord, user manual and test protocol) CURIS® 4 MHz Radiofrequenz-Generator (inkl. Netzkabel, Gebrauchsanweisung und Testprotokoll)
36 01 10	Foot switch two pedals for CURIS® (cut & coag), cable length: 4 m Zweipedal-Fußschalter CURIS®, Kabellänge: 4 m
37 01 54 L	Bipolar cable for CURIS®, length: 3 m Bipolares Kabel für CURIS®, Länge: 3 m
36 07 04	Monopolar handpiece (pencil) cut & coag, shaft 2.4 mm, cable length: 3 m Monopolarer Handgriff, Schaftdurchmesser 2,4 mm, Kabellänge: 3 m
36 02 38	Cable for single-use patient plates, length: 4.5 m Anschlusskabel für Einwegneutralelektroden, Länge: 4,5 m
36 02 22	Safety patient plates, single-use, packing: 5 x 10 pcs. Einwegneutralelektroden, Verpackung: 5 x 10 Stk.

Optional:

87 00 20	CURIS® basic set with reusable patient plate, cable length: 4 m (REF 36 02 26) CURIS® Basis-Set mit wiederverwendbarer Neutralelektrode, Kabellänge: 4 m (REF 36 02 26)
-----------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------



SUTTER MEDIZINTECHNIK GMBH

TULLASTRASSE 87 · 79108 FREIBURG/GERMANY · TEL. +49(0)761-51551-0 · FAX +49(0)761-51551-30