

MTM-1500 MED MTM-1400 MED

High-quality, low-maintenance treadmills



SCHILLER
The Art of Diagnostics

SCHILLER
MTM

MTM-1500 MED

SAFE – ACCURATE – QUIET


SCHILLER's treadmill ergometer MTM-1500 med is designed for exercise tests in sports medicine, cardiology, physiotherapy and fitness training. The ergometer is equipped and tested with the latest safety systems and has been manufactured to meet stringent quality requirements.

KEY FEATURES:

- ❖ Compatible with all SCHILLER exercise ECG devices
- ❖ With potential equalisation transformer for potential isolation from the mains, according to IEC 60601-1
- ❖ Automatic control by computer via RS-232 interface (interface cable 5 m included)
- ❖ Individually adjustable emergency stop button
- ❖ Steel tube hand rail; can be mounted on the right or left hand side. A second hand rail is optionally available.
- ❖ Automatic centring mechanism for the belt

ADVANTAGES:

- ❖ Very high durability thanks to outstanding manufacturing quality
- ❖ Excellent price / performance ratio
- ❖ Extremely safe
- ❖ Almost noiseless drive
- ❖ Maintenance-free motor
- ❖ Modern design
- ❖ Generation of reports with terminal software



MTM-1400 MED

THE NEW SCHILLER TREADMILL MTM-1400 MED

STABLE AND LOW-MAINTENANCE

The MTM-1400 med treadmill is virtually indestructible with its sturdy frame. It is very low maintenance and provides the runner or patient a comfortable run with its advanced construction. MTM-1400 med is distinguished by its smooth running, various functions, strong drive engine and the timeless design.

MEDICAL DEVICE TREADMILL

As a medical device, it is particularly suitable for use in the fields of cardiology, neurology, cardiac rehabilitation and physiotherapy.

The interface of the SCHILLER MTM-1400 med treadmill ergometer allows connection to ECG, ergospirometry systems, blood pressure measuring devices and software programs via the coscom® protocol.

CHILDREN HANDRAILS FOR PAEDIATRIC REHABILITATION

Children handrails are available as an additional option (only in combination with the safety arch). The MTM-1400 med treadmill can be used for paediatric rehabilitation with this option. Children can hold onto the handrail both on the sides and in front while walking on the treadmill, and therefore always have a sense of security and control.

NUMEROUS ADDITIONAL OPTIONS

You can customise the MTM-1400 med treadmill to your field of application – whether rehabilitation, cardiopulmonary diagnostics for patients or fitness – with the numerous additional options.





Americas
SCHILLER Americas Inc.
 Doral, Florida 33172

North America:
 Phone +1 786 845 06 20
 Fax +1 786 845 06 02
 sales@schilleramerica.com
 www.schilleramerica.com

Latin America & Caribbean
 Phone +1 305 591 11 21
 Fax +1 786 845 06 02
 sales@schillerla.ch
 www.schillerla.ch



Asia
SCHILLER Asia-Pacific / Malaysia
 52200 Kuala Lumpur, Malaysia
 Phone +603 6272 3033
 Fax +603 6272 2030
 sales@schiller.com.my
 www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
 A-4040 Linz
 Phone +43 732 70 99
 Fax +43 732 757 000
 sales@schiller.at
 www.schiller.at



China
Alfred Schiller (Beijing) Medical Equipment Co., Ltd.
 100015 Beijing, China
 Phone +86 010 52007020
 Fax +86 010 52007020-8016
 info@schillerchina.com
 www.schillermedical.cn



Croatia
SCHILLER d.o.o.
 10000 Zagreb
 Phone +385 1 309 66 59
 Fax +385 1 309 66 60
 info@schillerzg.hr
 www.schiller.hr



France
SCHILLER Médical S.A.S.
 F-67162 Wissembourg
 Phone +33 3 88 63 36 00
 Fax +33 3 88 63 36 49
 info@schiller.fr
 www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
 F-77608 Bussy St Georges
 Phone +33 1 64 66 50 00
 Fax +33 1 64 66 50 10
 contact@schillerfrance.fr
 www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
 D-85622 Feldkirchen b. München
 Phone +49 89 62 99 81 0
 Fax +49 89 609 50 90
 info@schillermed.de
 www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
 H-1141 Budapest
 Phone +36 1 383 4780 / +36 1 460 9491
 Fax +36 1 383 4778
 info@schillerhungary.hu
 www.schillerhungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
 Mumbai - 400 059, India
 Phone +91 22 6152 3333 / 2920 9141
 Fax +91 22 2920 9142
 sales@schillerindia.com
 www.schillerindia.com



Poland
SCHILLER Poland Sp. z o.o.
 PL-02-729 Warszawa
 Phone +48 22 843 20 89 / +48 22 647 35 90
 Fax +48 22 843 20 89
 schiller@schiller.pl
 www.schiller.pl



Serbia
SCHILLER d.o.o.
 11010 Beograd
 Phone +381 11 39 79 508
 Fax +381 11 39 79 518
 info@schiller.rs
 www.schiller.rs



Slovenia
SCHILLER d.o.o.
 2310 Slovenska Bistrica
 Phone +386 2 843 00 56
 Fax +386 2 843 00 57
 info@schiller.si
 www.schiller.si



Spain
SCHILLER ESPAÑA, S.A.
 ES-28232 Las Rozas/Madrid
 Phone +34 91 713 01 76
 Fax +34 91 355 79 33
 schiller@schiller.es
 www.schiller.es



Switzerland
SCHILLER-Reomed AG
 CH-8953 Dietikon
 Phone +41 44 744 30 00
 Fax +41 44 740 37 10
 sales@schiller-reomed.ch
 www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
 Okmeydanı-Sisli - Istanbul
 Phone +90 212 210 8681 (pbx)
 Fax +90 212 210 8684
 info@schiller.com.tr
 www.schiller-turkiye.com



United Kingdom
SCHILLER UK
 Bellshill, ML4 3PR
 Phone +44 1698 744 505
 Fax +44 1698 744 474
 sales@schilleruk.com
 www.schilleruk.com

